

Care Chair

Bariatric Stretcher / Recliner



CareChair 835

Integrated transfer chair, stretcher, and wheelchair. 31" seat width. Electric chair to stretcher conversion. Includes battery charger. 850lbs weight capacity.

Products Designed for Patient Safety in Transfer

The unique CareChair adjusts from full upright sitting position to full stretcher, and everything in-between. This allows the patient to gradually adjust from supine to upright, or to be put into a supine position for treatment or emergency care. The control is a simple battery powered switch on the back.

- Four models available to accommodate all patient sizes
- Central Locking Casters
- Infinite positioning adjustments, for improved pulmonary mechanics and patient comfort
- Foot and leg rests are separately adjustable
- Strategically located casters for superior stability, *Will not tip!*
- Easy storage with minimum footprint of 33" x 38"
- Adjustable height armrests
- Optional captive IV pole, oxygen tank holder, and integrated commode seat

The chair tilts back to redistribute the patient's weight, aiding in comfort and pressure relief. Chair also tilts forward to aid patient in egressing from chair. Stretcher has the ability to go into the trendelenberg position. Battery electric conversion.

Healthcare Rental Solutions

Specialty Medical
Equipment, LLC

Addressing The Needs of Every Patient

50 Kerry Place
Norwood, MA 02062
Tel. 781-844-5929
Fax. 781-297-8379

www.specialtymedicequipment.com

Care Chair

Bariatric Stretcher / Recliner



Model 835

Tilt Back



Specifications: CareChair Model 835

Weight Capacity	850 lbs
Seat Width	31 inches
Overall Width	39 inches

Seat Height Above Floor

Chair Position	24 inches
Stretcher Position	30 inches

Arm Rest (side rail) Height Above Seat

Chair Position 7-18 inches

Stretcher Position 1-10 inches

Recline Positioning Electrically Driven - Infinite

Battery 12 volt rechargeable sealed battery, good for approx. 75 cycles.
Recharge time is 12 hours

Foot Rests Adjustable individually, manual

Leg Rests Adjustable individually, manual and affected by other position changes

Options Oxygen Tank Holder,

Captive IV Rod and Holder (adjustable height)

Commode Seat

Shipping Weight Approximately 300 lbs (plus accessories)

Warranty 5 years workmanship and materials, 1 year on electrical system

Tilt Forward



Healthcare Rental Solutions

Specialty Medical
Equipment, LLC

Addressing The Needs of Every Patient

50 Kerry Place
Norwood, MA 02062
Tel. 781-844-5929
Fax. 781-297-8379

www.specialtymedicequipment.com