

Phoenix™ Knee CPM Series



The Phoenix™ Knee shown in upright easy storage position.

Designed With Hospital and In-Home Use In Mind

The Phoenix™ series of knee CPM is the result of McKelor™ Technologies' commitment to quality.

Reliable and versatile, the Phoenix™ is designed with in-home and hospital use in mind. Simple, easy-to-master hand held controls, combined with its lightweight design, the Phoenix™ is the perfect choice for your CPM needs.

Product Features

- Unique lift and bend motion provides a comfortable, accurate range of motion
- Controls are easily understood by staff and patients (1800's LED and 1850's dial controller)
- Units can be set up easily and quickly for small adult and pediatric patients without the need for additional kits
- May be programmed for flexion and/or extension pause
- Integrated bed mount system and non-skid pad prevent migration and maintain proper alignment
- Light weight, compact design make transportation and storage simple
- (Recommended) Phoenix™ Knee CPM Softgoods made exclusively by McKelor™ for the Phoenix™ Knee. Our softgoods feature comfortable, heavy pile and durable construction

Healthcare Rental Solutions

Specialty Medical
Equipment, LLC

Addressing The Needs of Every Patient

50 Kerry Place
Norwood, MA 02062
Tel. 781-844-5929
Fax. 781-297-8379

www.specialtymedicequipment.com

Phoenix™ Knee CPM Series



Specifications: Phoenix™ CPM Knee Series

Weight	26 lbs.
Length	28 Inches
Adjustable Length	42 Inches
Width	13 Inches
Range of Motion	-5° to 10°
Pause	0 to 30 Seconds
Limb Length	20 to 42 Inches
Calf Length	11 to 24 Inches
Thigh Length	9 to 18 Inches
Electrical Rating	14 V DC, 730 Amps
Warranty	Limited 1 year Parts & Labor



The **Phoenix™ Model 1850** offers a digital display pendant featuring patient lock-out, easy to use stop/start button and compliance meter. The hand-held control is easily understood and makes the hospital in-service and patient education simple.



The **Phoenix™ Model 1800** features an analog hand-held control pendant. The red control button allows the patient to stop or interrupt the unit as needed. Very simple and easy to use, this model is also ideal for in-home or hospital use.

Healthcare Rental Solutions

Specialty Medical
Equipment, LLC

Addressing The Needs of Every Patient

50 Kerry Place
Norwood, MA 02062
Tel. 781-844-5929
Fax. 781-297-8379

www.specialtymedicequipment.com