



Model Number: P6000

Dynamic Alternating Pressure Chair Pad Overlay

The Control Unit

Air Express™ A.P.M. Chair Unit

- Simple to use dial weight adjustment.
- Programmed never to exceed 32mmHg inflation pressure for optimal therapy.
- CPR quick hose release for rapid deflation.
- Compact and lightweight.
- Two year warranty

Air Express™ Economy A.P.M. Chair Unit

- Simple to use dial weight adjustment.
- Programmed never to exceed 32mmHg inflation pressure for optimal therapy.
- CPR quick hose release for rapid deflation.
- Compact and lightweight.
- One year warranty

The Therapeutic Support Surface

- Vapor permeable **FLWSHEAR II™ topcover**, stretchable and breathable to reduce maceration and friction. Machine washable.
- Patented differential cell design provides up to 4 inches of support in the crucial sacral region.
- Provides effective pressure management for patients using a geriatric chair.
- Overlay easily attaches to most models of geriatric chairs

Healthcare Rental Solutions

Specialty Medical
Equipment, LLC

Addressing The Needs of Every Patient

50 Kerry Place
Norwood, MA 02062
Tel. 781-844-5929
Fax. 781-297-8379

www.specialtymedicequipment.com

G.A.P Overlay



Model Number: P6000



Air Express™ A.P.M. P6000
Cost-effective system for dynamic alternating pressure relief.



The therapeutic support surface is lightweight and made of high quality materials. It conveniently rolls up for easy carrying with attached straps. Separate carrying case also available.

Specifications: G.A.P Geriatric Chair Overlay

Dimensions:	P6000D Air Express A.P.M. Control Unit: 5" x 10" x 12" / 12.7cm x 25.4cm x 30.5cm (5lbs. / 2.3kg) P6000I Air Express Economy A.P.M. Control Unit: 10" x 4" x 5" / 25.4cm x 10.2cm x 12.7cm (3.5lbs. / 1.6kg) Overlay: (2-4") x 19" x 70" / (5-10cm) x 48.3cm x 177.8cm (7lbs. / 3.2kg)
Agency Approvals:	P6000D: UL 2601, C22.2 No. 601-1-M90, TUV, EN 60601-1 P6000I UL544
AP Cycle Time:	P6000D Air Express A.P.M.: 5 Minutes P6000I Air Express Economy A.P.M.: 5 Minutes
Electrical Rating:	P6000D Air Express A.P.M.: 120V, 50Hz/60Hz, .15A P6000I Air Express Economy A.P.M.: 115V, 60Hz
Warranty	P6000D Air Express A.P.M. Control Unit: 2 years P6000I Air Express Economy A.P.M. Control Unit: 1 year Overlay: 1 year
Weight Limit	350lbs.

Healthcare Rental Solutions

Specialty Medical
Equipment, LLC

Addressing The Needs of Every Patient

50 Kerry Place
Norwood, MA 02062
Tel. 781-844-5929
Fax. 781-297-8379

www.specialtymedicequipment.com