


Liko FreeSpan Ultra Twin



Ultra Slingbar with safety connector.

Product Number:
303603 Ultra Slingbar



Ultra Twin on a FreeSpan railsystem - a freestanding overhead lift mounted on a multiple caster chassis. With a lifting capacity of 400kg (880 lbs.), Ultra Twin FreeSpan is designed for the heaviest lifting situations. The two interacting lift motors make it possible for the assistant to both lift and position the patient while in the sling during lift and transfer. The two lift motors may be simultaneously operated with one hand control, or used with separate hand controls.

FreeSpan Ultra Twin comes in two different models. For lifts with a sling or for horizontal lifts.

FreeSpan Ultra Twin is available with different rail lengths, ranging from 2.0-5.5 meters (6' 7½" - 18"). Choose length based on needs and the space available. FreeSpan UltraTwin is easy to assemble and move and the rails can be raised or lowered when needed.



Product Numbers

Products needed for a complete UltraTwin sling unit:

3103520	UltraTwin side supports	1 pair
3103115-50	FreeSpan Liftrail LR	2 pieces
3122007	Likorall 242 ES R2R	2 pieces
3126040	Ultra control unit	1 set
3126101-4	Battery charger (EU, UK US, AU)	1 piece
3103603	Ultra Slingbar	1 piece
	Appropriate sling/accessories	

Products needed for a complete UltraTwin stretcher unit:

3103521	UltraTwin sidesupport, wide	1 pair
3103115-50	FreeSpan Liftrail LR	2 pieces
3122007	Likorall 242 ES R2R	2 pieces
3126034	Handcontrol LR 242 ES, 2 button	2 pieces
3126101-4	Battery charger (EU, UK US, AU)	1 piece
3156058	Ultra Stretch	1 piece
	Appropriate sling/ accessories	

Maximum load: 400kg, 880lbs.

Healthcare Rental Solutions

Specialty Medical
Equipment, LLC

Addressing The Needs of Every Patient

50 Kerry Place
Norwood, MA 02062
Tel. 781-844-5929
Fax. 781-297-8379

www.specialtymedicequipment.com