



Titan-R/S  
Weight capacity 1000 lb.

## Standard Features

- Foot Egress bed that moves into a full chair position. Allowing you to enter and exit at the foot of the bed.
- Trapeze and inbed scale system included
- Fully electric DC motors, Rev. Trend, Trend, High/Low, Cardio Position
- Emergency battery back up
- 4 Padded drop down rails
- Therapeutic high density foam mattress
- Bed length 80"
- Bed deck low 19", high 34"
- Expandable Width - 36" deck - 54" deck

## Healthcare Rental Solutions

**Specialty Medical**  
Equipment, LLC

*Addressing The Needs of Every Patient*

50 Kerry Place  
Norwood, MA 02062  
Tel. 781-844-5929  
Fax. 781-297-8379

[www.specialtymedicequipment.com](http://www.specialtymedicequipment.com)