


Tri-Flex Bariatric Bed



Bariatric Bed
Shown with options



Optional Plexiglas Windowed
Head Rails


Multiple Feature High Performance Platform

Original Features of the Tri-Flex Bariatric Bed

- 1,000 lb. Capacity - Static and Transfer
- 2 Patient Surface Widths in One Frame
- 3 Piece Interchangeable Width Mattress
- "Fold & Roll" One Person Delivery System
- Integrated Scale System Option
- Exit/Enter at Foot of Bed Option

Second Generation Enhancements

- Sliding patient surface for easier expansion/retraction to 48" width
- Fold Down Foot Rails
- Plexiglas Windowed, Rotating Head Rails
- Narrower Frame for improved doorway clearance
- 1 Button Trendelenberg functions
- Improved Trapeze Efficiency with enhanced patient access
- New, *Streamlined* Dual Passage Assist
- Impulse Drive System



Compact "Fold & Roll"
1 Person Delivery System

Healthcare Rental Solutions

Specialty Medical
Equipment, LLC

Addressing The Needs of Every Patient

50 Kerry Place
Norwood, MA 02062
Tel. 781-844-5929
Fax. 781-297-8379

www.specialtymedicequipment.com

Tri-Flex Bariatric Bed



Improved Trapeze Option
for Greater Patient Access



3 & 4 Corner Setup Available



Angle Indicators

Easy to read and self setting; shows the degree of angle of the head gatch for recording levels of patient comfort relating to upper body position.



Push Indicators

Well positioned for optimum ease of use. Foam covered for a comfortable and safe grip. Affords less strain and effort for care givers while positioning or moving the bed. Also included with the Impulse Drive System.



Easy Patient Egress

Designed specifically for the Burke Tri-Flex II, this unique system enhancements provides the option of exiting or entering the bed through the foot end, encouraging greater patient mobility.



Burke Tri-Flex II Impulse Drive System

Safe and easy maneuverability through variable acceleration & positioning.

Healthcare Rental Solutions

Specialty Medical
Equipment, LLC

Addressing The Needs of Every Patient

50 Kerry Place
Norwood, MA 02062
Tel. 781-844-5929
Fax. 781-297-8379

www.specialtymedicequipment.com